

**PRODUCT BROCHURE**

# LAS-1000™

Liquid Autosampler for Seamless Integration  
into the XT Series Laboratory Analyzer  
for Total Sulfur-Total Nitrogen



## Precision

---

## Outstanding Repeatability

---

## Easy to Use

- Designed for ASTM D6667 and D4629 Applications
- Accurate and Precise Sample Handling
- XYZ-axis Vectorized Movement
- Fully Servo-motor Controlled
- Full Integration with XT Series Software
- Convenient Compact Footprint
- No Additional Benchtop Space Required

We offer the necessary tools and accessories to meet the growing needs and challenging demands of the analytical world. The LAS-1000™ Liquid Autosampler is built for ultra fast, reliable, efficient and safe handling of liquid samples and integrates seamlessly into the XT Series Laboratory Analyzer Total Sulfur – Total Nitrogen Solution.

The liquid autosampler does not require any additional benchtop space for operation. When combined with the XT Series it provides unmatched analytical performance making it highly versatile and suitable for a wide variety of applications and industries.

---

## Enhanced Features

The LAS-1000 is mounted directly on top of the XT Series analyzer eliminating the need for additional benchtop space. The use of servo-motors allows rapid, precise and accurate control of movements and eliminates routine maintenance needs. The use of XYZ-axis vectorized movement for sample handling reduces the time needed for operations.

- Large sample capacity: 162 Vials of 2 mL per tray holder with up to 2 trays holders
- Syringe protection with missing vial detection
- Optimized movement for highest sample analysis rate
- Highly reproducible injection
- Faster analysis for higher throughput

---

## Outstanding Repeatability

The LAS-1000, combined with the XT Series Analyzers, ensures optimal analytical performance. Autosampler injection repeatability below 0.60% RSD and linearity  $R > 0.9999$  allow the LAS-1000 combined with the XT Series to achieve 30 ppb repeatability for both total sulfur and total nitrogen measurement.

---


## Easy to Use

The LAS-1000 operation is fully integratable with the XT Series and managed directly through the Laboratory User Interface (LUI). The LUI software allows configuration and operation of the autosampler movements while giving users full control over the various analysis options.

\* On selected models

---

## 6 ppm Total Sulfur in Ultra-Low Sulfur Diesel



Run	1	2	3	4	5	6	7	8	9	10
Result (ppm)	6.03	6.04	6.07	5.99	6.04	5.99	6.03	5.99	6.00	6.02
Average 6.02			SD 0.025				RSD (%) 0.42			

## Outstanding Analytical Performance

The table above displays results of a 6.0 ppm/wt. Ultra Low Sulfur Diesel (ULSD) sample acquired with an XTS-1000 combined with the LAS-1000. The 10 run-graph above showcases the analytical performance and high speed of analysis of the XTS-1000/LAS-1000 system.

## SPECIFICATIONS

### Analytical Performance with XT-2000 Analyzer

**Measurement Method:** Total Sulfur – Excimer UV Fluorescence (EUVF)  
Total Nitrogen – Chemiluminescence

**ASTM Methods:** Total Sulfur – D5453  
Total Nitrogen – D4629

**Measurement Range:** Total Sulfur – 0-10%\*  
Total Nitrogen – 0-1%\*

**Repeatability:** Total Sulfur – 30 ppb or 2% whichever is greater  
Total Nitrogen – 30 ppb or 2% whichever is greater

**Analysis Time:** 3 minutes (typical)

**Typical Sample Size:** 5 to 20 µL

### Utility Requirements

**Power:** 100-240 VAC, 50/60 Hz

### Physical Parameters

**Dimensions:** 533 mm (21.0 in) Wide x 686 mm (27.0 in) Deep x 724 mm (28.5 in) High

**Weight:** 27 kg (60 lbs)

## GAIN REAL-TIME INSIGHT INTO YOUR PROCESS

Process Insights manufactures and delivers premium sensors, monitors, detectors, analyzers, instrumentation, and software that are mission-critical to keep your operations, personnel, and the environment safe – every day across the globe.

Get the most reliable, precision analytical technologies available on the market today. We will work to match your needs and budget, and provide the optimal, and most stable process analysis solution for your application.

---

## CENTERS OF EXCELLENCE | PROVIDING PROVEN SOLUTIONS

Process Insights is committed to solving our customers' most complex analytical, process, and measurement challenges everyday.

### Process Insights – The Americas

4140 World Houston Parkway Suite 180, Houston, TX 77032, USA +1 713 947 9591

### Process Insights – EMEA

ATRICOM, Lyoner Strasse 15, 60528 Frankfurt, Germany +49 69 20436910

### Process Insights – APAC

Wujiang Economic and Technology, Development Zone, No. 258 Yi He Road, 215200 Suzhou, Jiangsu Province, China +86 400 086 0106

---

For a complete range of products, applications, systems, and service options, please contact us at: [info@process-insights.com](mailto:info@process-insights.com)

For a complete list of sales & manufacturing sites, please visit: <https://www.process-insights.com/about-us/locations/>

COSA Xentaur, Tiger Optics, Extrel, Alpha Omega Instruments, ATOM Instrument, MBW Calibration, MGA, Guided Wave, ANALECT and LAR TOC Leader are trademarks of Process Insights, Inc.



REVOLUTIONIZING MEASUREMENT

**EVERYWHERE**